

Human CCL1 / C-C motif chemokine 1 ELISA Kit

Catalog No: #EK21037

Package Size: #EK21037 96T

Description

Product Name	Human CCL1 / C-C motif chemokine 1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Small-inducible cytokine A1,T lymphocyte-secreted protein I-309,SCYA1
Accession No.	Uniprot:P22362
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39pg/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.