

Human NPC1 / NPC intracellular cholesterol transporter 1 ELISA Kit

Catalog No: #EK21038

Package Size: #EK21038 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human NPC1 / NPC intracellular cholesterol transporter 1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Niemann-Pick C1 protein
Accession No.	Uniprot:O15118
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39pg/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.