

Human FAF1 / FAS-associated factor 1 ELISA Kit

Catalog No: #EK21065

Package Size: #EK21065 96T



Description

Product Name	Human FAF1 / FAS-associated factor 1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	hFAF1,UBX domain-containing protein 12,UBX domain-containing protein 3A,UBXD12,UBXN3A,CGI-03
Accession No.	Uniprot:Q9UNN5
Storage	6 months

Application Details

Detection Range:0.78-50ng/mLAssay Type:SandwichSensitivity:0.375ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.