

Human FKRP / Fukutin-related protein ELISA Kit

Catalog No: #EK21075

Package Size: #EK21075 96T



Description

Product Name	Human FKRP / Fukutin-related protein ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Ribitol-5-phosphate transferase
Accession No.	Uniprot:Q9H9S5
Storage	6 months

Application Details

Detection Range:0.312-20ng/mLAssay Type:SandwichSensitivity:0.13ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.