

Human TMEM201 / Transmembrane protein 201 ELISA Kit

Catalog No: #EK21089

Package Size: #EK21089 96T

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human TMEM201 / Transmembrane protein 201 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Spindle-associated membrane protein 1,NET5,SAMP1
Accession No.	Uniprot:Q5SNT2
Storage	6 months

Application Details

Detection Range:0.156-10ng/mLAssay Type:SandwichSensitivity:0.055ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.