

Human RPGR / X-linked retinitis pigmentosa GTPase regulator ELISA Kit

Catalog No: #EK21090

Package Size: #EK21090 96T

Description

Product Name	Human RPGR / X-linked retinitis pigmentosa GTPase regulator ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	RP3,XLRP3
Accession No.	Uniprot:Q92834
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.33ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.