

Human EPRS / Bifunctional glutamate/proline--tRNA ligase ELISA Kit

Catalog No: #EK21094

Package Size: #EK21094 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human EPRS / Bifunctional glutamate/proline--tRNA ligase ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Bifunctional aminoacyl-tRNA synthetase, Cell proliferation-inducing gene 32 protein, Glutamyl-prolyl-tRNA synthetase, GLNS, PARS, QARS, QPRS, PIG32
Accession No.	Uniprot:P07814
Storage	6 months

Application Details

Detection Range: 0.156-10ng/ml Assay Type: Sandwich Sensitivity: 0.055ng/ml

Note: This product is for in vitro research use only and is not intended for use in humans or animals.