

# Rat Csnk1a1 / Casein kinase I isoform alpha ELISA Kit

Catalog No: #EK21096

Package Size: #EK21096 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                       |
|--------------------|-------------------------------------------------------|
| Product Name       | Rat Csnk1a1 / Casein kinase I isoform alpha ELISA Kit |
| Applications       | ELISA                                                 |
| Species Reactivity | Rt                                                    |
| Other Names        | CKI-alpha,CK1                                         |
| Accession No.      | Uniprot:P97633                                        |
| Storage            | 6 months                                              |

## Application Details

Detection Range:78-5000pg/mlAssay Type:SandwichSensitivity:32pg/ml

Note: This product is for in vitro research use only and is not intended for use in humans or animals.