

Human DHRS3 / Short-chain dehydrogenase/reductase 3 ELISA Kit

Catalog No: #EK21097

Package Size: #EK21097 96T

Description

Product Name	Human DHRS3 / Short-chain dehydrogenase/reductase 3 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	DD83.1, Retinal short-chain dehydrogenase/reductase 1, retSDR1, Retinol dehydrogenase 17, Short chain dehydrogenase/reductase family 16C member 1, RDH17, SDR16C1, UNQ2424/PRO4983
Accession No.	Uniprot:O75911
Storage	6 months

Application Details

Detection Range: 0.312-20ng/ml Assay Type: Sandwich Sensitivity: 0.098ng/ml

Note: This product is for in vitro research use only and is not intended for use in humans or animals.