

Human BSG / Basigin ELISA Kit

Catalog No: #EK21099

Package Size: #EK21099 96T



Description

Product Name	Human BSG / Basigin ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	5F7, Collagenase stimulatory factor, Extracellular matrix metalloproteinase inducer, EMMPRIN, Leukocyte activation antigen M6, OK blood group antigen, Tumor cell-derived collagenase stimulatory factor, TCSF, CD147, UNQ6505/PRO21383
Accession No.	Uniprot:P35613
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:33pg/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.