

Human NPR1 / Atrial natriuretic peptide receptor 1
ELISA Kit

Catalog No: #EK21100

Package Size: #EK21100 96T



Description

Product Name	Human NPR1 / Atrial natriuretic peptide receptor 1 ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Atrial natriuretic peptide receptor type A,ANP-A,Guanylate cyclase A,GC-A,ANPRA
Accession No.	Uniprot:P16066
Storage	6 months

Application Details

Detection Range:0.312-20ng/mLAssay Type:SandwichSensitivity:0.15ng/mL

Note: This product is for in vitro research use only and is not intended for use in humans or animals.