

Mouse Sp1/ Transcription factor Sp1 ELISA Kit

Catalog No: #EK21105

Package Size: #319149 96T

Description

Product Name	Mouse Sp1/ Transcription factor Sp1 ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:78-5000pg/mL

Sensitivity:39.3pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.