

Human SYPL1/ Synaptophysin-like protein 1 ELISA Kit

Catalog No: #EK21109

Package Size: #EK21109 96T

Description

Product Name	Human SYPL1/ Synaptophysin-like protein 1 ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.78-50ng/mL

Sensitivity:0.345ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.