

Human BTN3A2/ Butyrophilin subfamily 3 member A2
ELISA Kit

Catalog No: #EK211111

Package Size: #EK211111 96T



Description

Product Name	Human BTN3A2/ Butyrophilin subfamily 3 member A2 ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78-5000pg/mL
Sensitivity:28pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.