

Mouse Pdgfra/ Platelet-derived growth factor receptor alpha ELISA Kit

Catalog No: #EK21124

Package Size: #EK21124 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Pdgfra/ Platelet-derived growth factor receptor alpha ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.171ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.