

Rat Eif2ak3/ Eukaryotic translation initiation factor 2-alpha kinase 3 ELISA Kit

Catalog No: #EK21135

Package Size: #EK21135-1 96T

Description

Product Name	Rat Eif2ak3/ Eukaryotic translation initiation factor 2-alpha kinase 3 ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Other Names	PRKR-like endoplasmic reticulum kinase, Pancreatic eIF2-alpha kinase, Pek, Perk

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.081ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.