

Human SPDL1/ Protein Spindly ELISA Kit

Catalog No: #EK21137

Package Size: #EK21137-1 96T

Description

Product Name	Human SPDL1/ Protein Spindly ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Arsenite-related gene 1 protein, Coiled-coil domain-containing protein 99, Rhabdomyosarcoma antigen MU-RMS-40.4A, Spindle apparatus coiled-coil domain-containing protein 1, hSpindly, CCDC99

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.083ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.