

Mouse Protein phosphatase 1 regulatory subunit 18(PPP1R18) ELISA Kit

Catalog No: #EK23340

Package Size: #EK23340-1 48T #EK23340-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Protein phosphatase 1 regulatory subunit 18(PPP1R18) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/mL Sensitivity:11.6pg/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.