

Human WD repeat-containing protein 1(WDR1) ELISA Kit

Catalog No: #EK23343

Package Size: #EK23343-1 48T #EK23343-2 96T

Description

Product Name	Human WD repeat-containing protein 1(WDR1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:27pg/mL

Sample Type:serum, plasma or other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.