

Mouse Pyruvate Dehydrogenase Kinase Isozyme 4 (PDK4) ELISA Kit

Catalog No: #EK23346

Package Size: #EK23346-1 48T #EK23346-2 96T



Description

Product Name	Mouse Pyruvate Dehydrogenase Kinase Isozyme 4 (PDK4) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.312-20ng/mL Sensitivity:0.112ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.