

Human MOB kinase activator 3C (MOB3C) ELISA Kit

Catalog No: #EK23347

Package Size: #EK23347-1 48T #EK23347-2 96T



Description

Product Name	Human MOB kinase activator 3C (MOB3C) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.105ng/mL Sensitivity:11.6pg/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.