

Human Serine protease HTRA3oO HTRA3oO ELISA Kit

Catalog No: #EK23349

Package Size: #EK23349-1 48T #EK23349-2 96T

Description

Product Name	Human Serine protease HTRA3oO HTRA3oO ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL Sensitivity:0.062ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.