

Human UTP--glucose-1-phosphate uridylyltransferase (UGP2) ELISA Kit

Catalog No: #EK23350

Package Size: #EK23350-1 48T #EK23350-2 96T

Description

Product Name	Human UTP--glucose-1-phosphate uridylyltransferase (UGP2) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL Sensitivity:0.057ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50 μ L Assay Time:1-4.5h
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.