

# Human PAS domain-containing serine/threonine-protein kinase(PASK) ELISA Kit

Catalog No: #EK23353

Package Size: #EK23353-1 48T #EK23353-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                             |
|--------------------|-----------------------------------------------------------------------------|
| Product Name       | Human PAS domain-containing serine/threonine-protein kinase(PASK) ELISA Kit |
| Applications       | ELISA                                                                       |
| Species Reactivity | Human                                                                       |

## Application Details

Detect Range:0.781-50ng/mL Sensitivity:0.32ng/mL Sample Type:serum, plasma or other biological fluids Sample Volume: 50  $\mu$ L Assay Time:1-4.5h  
Detection wavelength:450 nmDetection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.