

APC anti-human HLA-A*01:01 B*08:01 C*08:01 HLA Class II

Catalog No: #FC0112

Package Size: #FC0112 25 Tests

Description


Product Name	APC anti-human HLA-A*01:01 B*08:01 C*08:01 HLA Class II
Clonality	Monoclonal
Clone No.	s-W6/32
Isotype	mouse IgG2a, κ
Applications	FC
Species Reactivity	human
Conjugates	APC
Target Name	HLA-A*01:01 B*08:01 C*08:01 HLA Class II
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.


Application Details

Test: 5 μl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images





Human peripheral blood lymphocytes stained with APC anti-human HLA Class I (colored histogram) or APC mouse IgG2a κ (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.