

## PE-Cy5 anti-human/mouse CD11b

Catalog No: #FC0161

Package Size: #FC0161-1 25Tests #FC0161-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE-Cy5 anti-human/mouse CD11b                                     |
| Host Species       | rat                                                               |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-M1/70                                                           |
| Isotype            | rat IgG2b                                                         |
| Applications       | FC                                                                |
| Species Reactivity | human/mouse                                                       |
| Conjugates         | PE-Cy5                                                            |
| Target Name        | CD11b                                                             |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood granulocytes stained with PE-Cy5 anti-human/mouse CD11b (colored histogram) or PE-Cy5rat IgG2b (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.