

PE-Cy5 anti-human CD14

Catalog No: #FC0162

Package Size: #FC0162-1 25Tests #FC0162-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PE-Cy5 anti-human CD14
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-UCHM-1
Isotype	mouse IgG2a
Applications	FC
Species Reactivity	human
Conjugates	PE-Cy5
Target Name	CD14
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood monocytes stained with PE-Cy5 anti-human CD14 (colored histogram) or PE-Cy5 mouse IgG2a (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.