

FITC anti-human CD154

Catalog No: #FC0314

Package Size: #FC0314-1 25Tests #FC0314-2 100Tests

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description


Product Name	FITC anti-human CD154
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-M91
Isotype	mouse IgG1
Applications	FC
Species Reactivity	human
Conjugates	FITC
Target Name	CD154
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Intracellular Cell Stimulation Cocktail -stimulated human PBMCs (4hours) stained with FITC anti-human CD154 (colored histogram) or FITC mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.