

FITC anti-human HLA-AoO BoO CoO HLA Class loO

Catalog No: #FC0328

Package Size: #FC0328 25 Tests

Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | FITC anti-human HLA-AoO BoO CoO HLA Class loO                     |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-W6/32                                                           |
| Isotype            | mouse IgG2a, κ                                                    |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | FITC                                                              |
| Target Name        | HLA-AoO½oO½BoO½oO½CoO½oO½HLA Class loO½oO½                        |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |


Application Details

Test: 5 &mu;l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images





Human peripheral blood lymphocytes stained with APC anti-human HLA Class I (colored histogram) or APC mouse IgG2a  $\kappa$  (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.