

FITC anti-human HLA-DQ (MHC Class II)

Catalog No: #FC0331

Package Size: #FC0331-1 25Tests #FC0331-2 100Tests

Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | FITC anti-human HLA-DQ (MHC Class II)                             |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-1a3                                                             |
| Isotype            | mouse IgG2a                                                       |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | FITC                                                              |
| Target Name        | HLA-DQ (MHC Class II)                                             |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with FITC anti-human HLA-DQ (colored histogram) or FITC mouse IgG2a (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.