

FITC anti-HLA-G

Catalog No: #FC0335

Package Size: #FC0335-1 25Tests #FC0335-2 100Tests

Description

Product Name	FITC anti-HLA-G
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-87G
Isotype	mouse IgG2a, κ
Applications	FC
Conjugates	FITC
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1×10^6 cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.