

FITC anti-human Ki-67

Catalog No: #FC0343

Package Size: #FC0343-1 25Tests #FC0343-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	FITC anti-human Ki-67
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-3D5
Isotype	mouse IgG2a
Applications	FC
Species Reactivity	human
Conjugates	FITC
Target Name	Ki-67
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



MCF-7 cell line was fixed and permeabilized stained with FITC anti-human Ki-67 Antibody (colored histogram) or FITC mouse IgG2a (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.