

## Alexa Fluor 488 anti-human CD3

Catalog No: #FC0351

Package Size: #FC0351-1 25Tests #FC0351-2 100Tests

## Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | Alexa Fluor 488 anti-human CD3                                    |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-OKT3                                                            |
| Isotype            | mouse IgG2a                                                       |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | AF Dye                                                            |
| Target Name        | CD3                                                               |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes stained with Alexa Fluor 488 anti human CD3 (colored histogram) or Alexa Fluor 488 mouse IgG2a (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.