

## Alexa Fluor 488 anti-human CD16

Catalog No: #FC0355

Package Size: #FC0355-1 25Tests #FC0355-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | Alexa Fluor 488 anti-human CD16                                   |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-3G8                                                             |
| Isotype            | mouse IgG1,k                                                      |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | AF Dye                                                            |
| Target Name        | CD16                                                              |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes stained with Alexa Fluor 488 anti-human CD16 (colored histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.