

PE-Cy7 anti-human CD274

Catalog No: #FC0429

Package Size: #FC0429-1 25Tests #FC0429-2 100Tests



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PE-Cy7 anti-human CD274
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-29E2A3
Isotype	mouse IgG2b, κ
Applications	FC
Species Reactivity	human
Conjugates	PE-Cy7
Target Name	CD274
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



PHA-stimulated (3 days) human peripheral blood lymphocytesstained with PE-Cy7 anti-human CD274 (colored histogram) or PE-Cy7 mouse IgG2b κ (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.