

## PE anti-human CD10

Catalog No: #FC0475

Package Size: #FC0475-1 25Tests #FC0475-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE anti-human CD10                                                |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-CB-CALLA                                                        |
| Isotype            | mouse IgG2b, $\kappa$                                             |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | PE                                                                |
| Target Name        | CD10                                                              |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood granulocytes stained with PE anti-human CD10 (colored histogram) or PE mouse IgG2b  $\kappa$  (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.