

PE anti-human CD49f

Catalog No: #FC0535

Package Size: #FC0535-1 25Tests #FC0535-2 100Tests



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PE anti-human CD49f
Host Species	rat
Clonality	Monoclonal
Clone No.	s-4F10
Isotype	rat IgG2b, κ
Applications	FC
Species Reactivity	human
Conjugates	PE
Target Name	CD49f
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with PE anti-human CD49f (colored histogram) or PE mouse IgG2b κ (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.