

PE anti-human CD51

Catalog No: #FC0537

Package Size: #FC0537-1 25Tests #FC0537-2 100Tests



**SAB**<sup>®</sup>

Signalway Antibody

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE anti-human CD51                                                |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-13C2                                                            |
| Isotype            | mouse IgG1                                                        |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | PE                                                                |
| Target Name        | CD51                                                              |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



MCF-7 cell line stained with PE anti-human CD51 (colored histogram) or PE mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.