

PE anti-human CD127

Catalog No: #FC0586

Package Size: #FC0586-1 25Tests #FC0586-2 100Tests

Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE anti-human CD127                                               |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-A019D5                                                          |
| Isotype            | mouse IgG1,k                                                      |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | PE                                                                |
| Target Name        | CD127                                                             |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood CD4+ T cells stained with PE anti-human CD127 and APC anti-human CD25

Note: This product is for in vitro research use only and is not intended for use in humans or animals.