

## PE anti-human CD307b

Catalog No: #FC0613

Package Size: #FC0613-1 25Tests #FC0613-2 100Tests

## Description

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE anti-human CD307b                                              |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-b24                                                             |
| Isotype            | mouse IgG2a, $\kappa$                                             |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | PE                                                                |
| Target Name        | CD307b                                                            |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.