

PE anti-human HLA-DP (MHC Class II) Monomorphic

Catalog No: #FC0623

Package Size: #FC0623-1 25Tests #FC0623-2 100Tests

Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE anti-human HLA-DP (MHC Class II) Monomorphic                   |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-B7/21                                                           |
| Isotype            | mouse IgG3                                                        |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | PE                                                                |
| Target Name        | HLA-DP (MHC Class II) Monomorphic                                 |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with PE anti-human HLA-DP (colored histogram) or PE mouse IgG3 (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.