

PE-Cy5.5 anti-human CD3

Catalog No: #FC0642

Package Size: #FC0642-1 25Tests #FC0642-2 100Tests

Description


Product Name	PE-Cy5.5 anti-human CD3
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-UCHT1
Isotype	mouse IgG1,k
Applications	FC
Species Reactivity	human
Conjugates	PE-Cy5.5
Target Name	CD3
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with PE-CY5.5 anti-human CD3 (colored histogram) or PE-CY5.5 mouse IgG1 κ (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.