

PerCP-Cy5.5 anti-human CD8

Catalog No: #FC0678

Package Size: #FC0678-1 25Tests #FC0678-2 100Tests



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PerCP-Cy5.5 anti-human CD8
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-MEM-31
Isotype	mouse IgG2a
Applications	FC
Species Reactivity	human
Conjugates	PerCP-Cy5.5
Target Name	CD8
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with PerCP-Cy5.5 anti-human CD8 Antibody (colored histogram) or PerCP-Cy5.5 mouse IgG2a (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.