

Pacific Blue anti-human CD65

Catalog No: #FC0737

Package Size: #FC0737-1 25Tests #FC0737-2 100Tests

Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | Pacific Blue anti-human CD65                                      |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-VIM8                                                            |
| Isotype            | mouse IgM, κ                                                      |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | Pacific Blue                                                      |
| Target Name        | CD65                                                              |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood stained with Pacific Blue anti-human CD65

Note: This product is for in vitro research use only and is not intended for use in humans or animals.