

Pacific Blue anti-human CD18

Catalog No: #FC0788

Package Size: #FC0788-1 25Tests #FC0788-2 100Tests

Description


Product Name	Pacific Blue anti-human CD18
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-TS1/18.2.11
Isotype	mouse IgG1 $\alpha\Omega\%o\Omega\%k$
Applications	FC
Species Reactivity	human
Conjugates	Pacific Blue
Target Name	CD18
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with Pacific Blue anti-human CD18 (colored histogram) or Pacific Blue mouse IgG1 κ (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.