

## PE anti-human FOXP3

Catalog No: #FC0801

Package Size: #FC0801-1 25Tests #FC0801-2 100Tests

## Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE anti-human FOXP3                                               |
| Host Species       | rat                                                               |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-PCH101                                                          |
| Isotype            | rat IgG2a, k                                                      |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | PE                                                                |
| Target Name        | FOXP3                                                             |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes was fixed permeabilized and intracellular stained with PE anti-human FOXP3 and APC anti-human CD25

Note: This product is for in vitro research use only and is not intended for use in humans or animals.