

## PE-Cy5.5 anti-human /mouse CD45R/B220

Catalog No: #FC1192

Package Size: #FC1192-1 25Tests #FC1192-2 100Tests

## Description

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | PE-Cy5.5 anti-human /mouse CD45R/B220                             |
| Host Species       | rat                                                               |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-RA3-6B2                                                         |
| Isotype            | rat IgG2a                                                         |
| Applications       | FC                                                                |
| Species Reactivity | human /mouse                                                      |
| Conjugates         | PE-Cy5.5                                                          |
| Target Name        | CD45R/B220                                                        |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.