

Purified Armenian Hamster IgG isotype Control

Catalog No: #FC1235

Package Size: #FC1235-1 100ug #FC1235-2 1mg



Description

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | Purified Armenian Hamster IgG isotype Control                     |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-PIP                                                             |
| Applications       | FC                                                                |
| Species Reactivity | Armenian Hamster                                                  |
| Conjugates         | Purified                                                          |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.