

# FITC Mouse IgG2a, $\kappa$ isotype Control

Catalog No: #FC1247

Package Size: #FC1247-1 25Tests #FC1247-2 100Tests

## Description

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | FITC Mouse IgG2a, $\kappa$ isotype Control                        |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-MOPC-173                                                        |
| Applications       | FC                                                                |
| Species Reactivity | Mouse                                                             |
| Conjugates         | FITC                                                              |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5  $\mu$ l/TestNotice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically,  $1 \times 10^6$  cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.