

Pacific Blue anti-mouse IgM Isotype Control

Catalog No: #FC1297

Package Size: #FC1297-1 25Tests #FC1297-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pacific Blue anti-mouse IgM Isotype Control
Clonality	Monoclonal
Clone No.	s-7C2C6C12
Applications	FC
Species Reactivity	Mouse
Conjugates	Pacific Blue
Target Name	IgM Isotype Control
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.